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FACSAria III is a cytometer analyzer separator of cells suspensions of high speed and resolution, capable of analyzing up to 13 parameters at the same time, 11 of them fluorescent (5B + 3R + • 3V).

The automatic adjustment of the drop delay and the automatic monitoring of the software classification, allows to establish a delay in the fall and thus ensure the obtaining of a high reproduction in the separation, and the remarkable improvement of the purity.

#### **Excitation and detection:**

##### **488 nm laser:**

- 530/30-nm FITC, GFP
- 582/15-nm PE
- 610/20-nm ,PE-Texas Red<sup>®</sup>,
- 695/40-nm PerCP-Cy5.5 ,
- 780/60-nm PE-Cy7

##### **633 nm laser:**

- APC of 660/20 nm
- Alexa fluor<sup>®</sup> of 730/45-nm
- 780/60-nm APC-Cy7 or APC-H7

##### **405 nm laser:**

- Blu Ray Brilliant 421 BD 450/50-nm, BD Horizon <sup>™</sup> V450, Pacific Blue <sup>™</sup>, DAPI
- 510/50-nm BD Brilliant Violet 510, BD Horizon V500, AmCyan
- 610/20-nm BD Brilliant Violet 605

Optical detection system of octagons and trigons based on the initial detection of long wavelengths, less energetic thus avoiding signal loss

**Sensitivity:** FITC: 85 MESF, PE: 29 MESF • Sensitivity in FSC and SSC: separation of 0.5 micron microspheres of noise in lateral dispersion.

**Resolution** in front light scattering is optimized for the detection of lymphocytes, monocytes and granulocytes

#### **Electronic System**

It is a fully digital system of all signals eliminating electronic dead time, thus increasing the real cellular performance. Total intra and interlaser compensation for all parameters, both "online" and "offline".

#### **Cell Sorting System**

Four-way sorting, cooling / heating system for tubes or sample collection plates, cloning unit (ACDU) for automatic separation of cells in plates slides or in eppendorf tubes, sample collection area separated and protected. It has nozzle of 70, 80 and 100

microns for the different sizes of particles. It allows the separation at high speed with purities greater than 98% and maintaining viability and cellular functionality.